


The efficacy of  
a combined  
action



**Imperium<sup>®</sup>**  
**master**



with  
**RF FRACTIONAL**  
NEW TOOL



# Imperium<sup>®</sup> master

**Imperium Master** is dedicated to the aesthetic treatment of the face. The approach may be Endogenous (Biostimulation) and Exogenous (Transdermal delivery of formulas).

**Imperium Master** is not just an ordinary beauty treatment: it is an innovative technology that helps your body regenerate itself in a natural manner, immediately improving your aspect and the overall well-being.


With **Imperium Master** results are immediately seen: your face looks young and bright.

## FIELDS OF APPLICATION

- Microroughness
- Acne scars
- Medium, moderate Wrinkles
- Stretch Marks
- Enlarged pores
- Periorbital wrinkles and Ptosis

## IMPERIUM MASTER ADOPTS THE EXCLUSIVE TECHNIQUE WITH RETICULAR FRACTIONAL FPC TIP

The TIP creates 125 micro holes in the epidermis, leading to a simultaneous denaturation of collagen in the dermis due to the emission of pulsed radiofrequency. The process involves the deep layers of the dermis, preserving the superficial ones. The tissues are regenerated in a completely natural manner, with no downtime.




Before

FPC

After


The thermal stimulation of the deep dermis, together with mechanical insult of the perforation, leads to a strong collagen remodeling and a whole series of benefits and advantages:



steel tips  
(AISI 316L)

**FPC Papillary  
Fractional  
Pulsed  
Crown**

- Fibroblast Skin Tightening
  - True skin rejuvenation
  - Reduction of deep wrinkles
  - Development of new collagen
  - Firming and remodeling of the dermis
- 
- Painless treatment
  - Quick recovery time after the session
  - Absence of an effect on melanin which makes the procedure safe in darker skin types and on lighter targets
  - Effective already since the first two sessions.
  - Few sessions (4-5) spaced by 3/4 weeks for rapid and lasting results



## SONOBLADE

**EPIDERMOABRASION** FOR AN EFFICACIOUS INTENSE PEELING. **SONOPHORESIS** FOR ABSORPTIONS OF SERUMS

The corneous layer is a protective barrier for the body and it acts as a filter to prevent penetration of numerous substance, whatever is their molecular weight.

**Imperium Master** has a function that removes up to 50% of the protective corneous layer in a non invasive way. This treatment is useful for increasing cutaneous permeability and making more fluid the hydrolypidic film.

Its ultrasonic action at 25 Khz also allows the transdermal delivery of serums and specific fluid formulations.

## RF SPHERES

- Biostimulation of the face, neck, décolleté
- Biostimulation of fibroblast activity, collagen (III) and Elastin
- Biostimulation of hyaluronic acid
- Oxygenation of tissues
- Deep vascularization


The controlled heating of the derma causes the contraction of collagen fibers, while the action inflicted on tissues by the heat stimulates the proliferation of Fibroblasts, thereby activating endogenous tissue regeneration.

This makes **Imperium MASTER** a unique system broadening its range of applications and significantly reducing regeneration time compared to the conventional equipments, facilitating by far the mode of application thanks also to a series of actuator of advanced technology.



## RF MULTIPOLAR

For the treatments of **wide areas** with an integrated system of **vibropercussion** and temperature control.



## VIBROPERCUSSION

**Imperium Master** system has special integrated and **patented actuators** which combine the emission of physical energies with a **percussion vibration system**, an ideal synergy for deep and highly stimulating endothermal massage.

This **exclusive characteristic** increases tolerance to heat, blood flow and lymphatic drainage.

**CONTACT US NOW**  
**0974 19.36.516**

[www.breramedical.com](http://www.breramedical.com)  
[info@breramedical.com](mailto:info@breramedical.com)

DISTRIBUTOR

**brera**<sup>®</sup>  
medical technologies

COMPANY WITH CERTIFIED QUALITY MANAGEMENT SYSTEM  
**UNI EN ISO 9001:2015**  
Patented n. 0000281819 - 0000275485 - 0000270042 - 0000281819

The data, images and all other described here, have no contract value  
and may be modified at any time without prior notice.

Imperium<sup>®</sup> and Brera<sup>®</sup>  
are registered trademarks and may not be produced.

MADE IN ITALY

[www.breramedical.com](http://www.breramedical.com)