

Imperium®



the efficacy of a combined action

Imperium Master is dedicated to the aesthetic treatment of the face. The therapeutic approach may be Endogenous (Biostimulation) and Exogenous (Transdermal delivery of formulas without using needles).

Imperium Master is not just an ordinary beauty treatment: it is an innovative technology that helps your body regenerate itself in a natural manner, immediately improving your aspect and the overall well-being. With **Imperium Master** results are immediately seen: your face looks young and bright.

2018 New combination: Imperium Master adopt the unique **BRERA reticular fractional technique with crown tip**, a bipolar configuration with three-dimensional effect.

The applicator creates 225 micro holes in the epidermis, leading to a simultaneous denaturation of collagen in the dermis due to the emission of pulsed radiofrequency.

FIELDS OF APPLICATION

- Microroughness
- Acne scars
- Medium, moderate Wrinkles
- Stretch Marks
- Enlarged pores
- Periorbital wrinkles and Ptosis

HANDPIECES




FPC Papillary RF

Innovative fractionated handpiece, with 225 micro-needles on a square area of 1.2 x 1.2 cm. The handpiece is designed to perform the treatment with tips at controlled pressure (through an exclusive cushioning system). It allows to create micro holes at various depths in a uniform way, reducing pain and eliminating any chance of bleeding.



Sonoblade

Epidermoabrasion.
Sonophoresis.
For an efficacious intense peeling.




RF Sphere

Bipolar treatment for jaws, wrinkles, forehead and hands.




RF Multipolar

For the treatments of wide areas with an integrated system of vibropercussion and temperature control.



FPC Papillary Fractional Pulsed Crown



Before Treatment Imperium MASTER FPC Papillary Treatment After Treatment

The thermal stimulation of the deep dermis, together with mechanical insult of the perforation, leads to a strong collagen remodeling and a whole series of benefits and advantages:

- True skin rejuvenation and tightening
- Reduction of deep wrinkles
- Development of new collagen
- Firming and remodelling of the dermis
- Painless treatment, without anesthesia
- Quick recovery time after the session
- Absence of an effect on melanin which makes the procedure safe on darker skin types
- Effective already since the first two sessions
- Few sessions (4-5) spaced by 3/4 weeks for rapid and lasting results

brera[®]
medical technologies

CONTACT US NOW
+390974 19.36.516

UNI CEI EN ISO 13485:2012 - UNI EN ISO 9001:2008 - Patented n. 0000275485

The data, images and all other described here, have no contract value and may be modified at any time without prior notice.

www.breramedical.com

CE
1370

The FPC Papillary handpiece acts simultaneously on the epidermis and dermis without side effects and without long recovery times after the treatment.

The efficacy of
a combined
action



Imperium[®]
master



with
RF FRACTIONAL
NEW TOOL

brera[®]
medical technologies